



BioStep Inc. 7221 Clybourn Ave.  
Sun Valley, CA 91352  
818 373 0010 T  
818 373 0030 F  
www.biosteportho.com

# FUNCTIONAL FOOT ORTHOTIC RX FORM

DATE: \_\_\_\_\_

## CLIENT INFORMATION

## BILLING INFO

## SHIPPING INFO

SAME AS BILLING INFO?

NAME \_\_\_\_\_  
ADDRESS \_\_\_\_\_  
CITY \_\_\_\_\_ STATE \_\_\_\_\_ ZIPCODE \_\_\_\_\_  
PHONE (\_\_\_\_) \_\_\_\_ - \_\_\_\_\_ FAX (\_\_\_\_) \_\_\_\_ - \_\_\_\_\_  
EMAIL \_\_\_\_\_ PO# \_\_\_\_\_  
CASTING CONTACT \_\_\_\_\_

NAME \_\_\_\_\_  
ADDRESS \_\_\_\_\_  
CITY \_\_\_\_\_ STATE \_\_\_\_\_ ZIPCODE \_\_\_\_\_  
PHONE (\_\_\_\_) \_\_\_\_ - \_\_\_\_\_ FAX (\_\_\_\_) \_\_\_\_ - \_\_\_\_\_

STANDARD DELIVERY

RUSH DELIVERY + EXTRA CHARGE

Needed By: \_\_\_\_/\_\_\_\_/\_\_\_\_

## PATIENT INFO

NAME \_\_\_\_\_

GENDER \_\_\_\_\_ WEIGHT \_\_\_\_\_ AGE \_\_\_\_\_ HEIGHT (Inches) \_\_\_\_\_ \*SHOE SIZE (Required) \_\_\_\_\_

SHOE STYLE \_\_\_\_\_  SHOES ENCLOSED  INSOLES ENCLOSED  TRACING ENCLOSED

## STEP 1

### CHOOSE A BASE ORTHOTIC TYPE

#### FUNCTIONAL ORTHOTICS

- Rigid - 4mm Poly, Rigid
- Multipurpose - 3mm Poly, Semirigid

#### SPORT ORTHOTICS

- Running - 3mm Poly, 30 Varus FF Running Wedge
- Aerobic - 2mm Poly, Soft EVA Arch fill, Flexible
- Basketball - 3mm Poly, Soft EVA Arch fill, Semirigid
- Sport Graphite - Carbon shell, Semirigid

#### ACCOMMODATIVE ORTHOTICS

- Geriatric - 2mm Poly, Semi-Flexible
- Bio-Cork - Thermocork Shell, Soft Covers Standard.
- BioStep TLCO - Deep heel cup, high Medial/Lateral walls.  
All Cork and EVA Construction

Diabetic Insert - See our "Custom Diabetic Order Form"

#### PEDIATRIC ORTHOTICS

- Gait Plate - Lateral Wall, Deep Heel Cup  
Please Specify: Promote In-Toeing or Out-Toeing
- U.C.B.L. - 3mm Poly, Deep Heel Cup with lateral wall medial flare.

#### FASHION ORTHOTICS

- Fashion Graphite - Flat Heel Cup, Graphite Shell, No Post
- Poly-Fashion - Flat Heel Cup, Poly Shell, No Post

## STEP 2

### CHOOSE SHELL ADDITIONS

L  R

- Morton's Extension
- Reverse Morton's Extension
- Medial Flare
- Lateral Wall
- Heel Punch

## STEP 3

### CHOOSE POSTING

YES EXTRINSIC  NO POST  NEUTRAL INTRINSIC

#### Extrinsic Rear Foot

Varus  L \_\_\_\_\_  R \_\_\_\_\_  
Valgus  L \_\_\_\_\_  R \_\_\_\_\_

#### Extrinsic Fore Foot

Varus  L \_\_\_\_\_  R \_\_\_\_\_  
Valgus  L \_\_\_\_\_  R \_\_\_\_\_

#### Intrinsic Rear Foot

Varus  L \_\_\_\_\_  R \_\_\_\_\_  
Valgus  L \_\_\_\_\_  R \_\_\_\_\_  
Heel Lift  L \_\_\_\_\_ mm

#### Intrinsic Fore Foot

Varus  L \_\_\_\_\_  R \_\_\_\_\_  
Valgus  L \_\_\_\_\_  R \_\_\_\_\_  
mm  R \_\_\_\_\_ mm

## STEP 4

### CHOOSE PADS AND ACCOMODATIONS

L  R

- Met Pad 2-4
- Met Bar 1-5
- Heel Pad
- Horse Shoe Pad
- Dancers Pad

MPJ Accommodations					
	1	2	3	4	5
L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
R	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>



CONTINUED ON BACK OF THIS PAGE



**STEP 5**

**CHOOSE COVERS**

**OPTIONAL COVERS (PLUS CHARGE)**

**Length**

- Mets
- Sulcus
- Toes

**Standard Covers**

- Simulated Leather
- Spenco
- Plastazote
- Perforated EVA
- Non-Perforated EVA
- Swirl EVA (Choose Color)
  - Smoke Swirl
  - Green Swirl
  - Blue Swirl

**Genuine Leather**

- Black
- Brown

**Perforated Genuine Leather**

- Midnight Black
- Steel Gray
- Toffee
- Moonlight Blue
- Angel Blue
- Rosewater

**Total Cover Thickness:**

1/16" 1/8" 3/16" 1/4"

**STEP 6**

**SPECIAL INSTRUCTIONS**

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---